

■ Specification

Model	Processing flow rate (L/h at 10°C)	Dimension (mm)			Operating weight (kg)	Power supply (V)	Power consumption (VA)	Pure water tank (L)	Quality of treated water
		W	D	H					
PR-0040-SG	40 (0.7L/min)	350	639	1142	98	AC100 ± 10%	450	20	Below 1 μ S/cm
PR-0100-SG	100 (1.7L/min)				100		500		
PR-0250-SG	250 (4.2L/min)	520	908	1400	220		900	100	
PR-0500-SG	500 (8.3L/min)				250		1300		