

■ Specification

| Model | Processing flow rate (L/h at25°C) | Dimension (mm) | | | Operating weight (kg) | Power supply (V) | Power capacity (VA) | Pure water tank (L) | Quality of treated water |
|------------|--------------------------------------|----------------|-----|------|--------------------------|------------------------|------------------------|------------------------|--------------------------|
| | | W | D | H | | | | | |
| SD-0040-SG | 40 | 400 | 562 | 1500 | Approx.110 | AC100±10% (50/60Hz) | 600 | 20 | Above 10MΩ·cm |
| SD-0080-SG | 80 | 450 | 804 | 1480 | Approx.175 | | 1500 | 60 | |
| SD-0200-SG | 200 | 600 | 951 | 1500 | Approx.280 | | | 100 | |