

■ Specification

Model		Processing flow rate (L/h)	Dimension (mm)			Operating weight (kg)	Power supply (V)	Power capacity (kVA)	Final filter	UV	TOC monitor	Quality of treated water
			W	D	H							
FP-0500	SM-000	500	350	750	1400	Approx.100	AC100±10% 50-60Hz universal	1	MF	Sterilization	-	Above 18.0MΩ·cm
	SU-000	500						1	UF	Sterilization	-	
	XM-0T0	500						1	MF	Oxidation	○	
	XU-0T0	500						1	UF	Oxidation	○	