

■ Specification

Model	Quality of treated water	TOC (ppb)	Processing flow rate	Dimension (mm)			Operating weight (kg)	Power consumption (VA)	Final filter	TOC monitor
				W	D	H				
FP-0120 $\alpha$ -UT1	18.2 M $\Omega$ · cm	1-3	< 2L/min	354	335	448	Approx.26	130 (AC100-240V)	UF	Real time
FP-0120 $\alpha$ -MT1									MF	
FP-0120 $\alpha$ -M01		-								-