

■ Specification

Model	Quality of treated water	Processing flow rate	Dimension (mm)			Operating weight	Power consumption	Final filter	Purewater Tank
			W	D	H	(kg)	(VA)		
UP-0090 α -TU1	18.2 M Ω ·cm	< 1.5L/min	354	446	460	Approx.32	130 (AC100-240V)	UF	Built In
UP-0090 α -TM1								MF	Built In
UP-0090 α -0U1			354	335	448	Approx.28		UF	Option
UP-0090 α -0M1								MF	Option

■ Option Purewater Tank(for 0U1/0M1)

Tank Volume (L)	Dimension (mm)			Operating weight
	W	D	H	(kg)
20	255	455	542	約32
30	465	419	826	約47
60	465	419	1091	約77
100	465	419	1526	約122