

■ Specification

Model	Quality of treated water※	TOC※	Processing flow rate	Dimension (mm)			Operating weight (kg)	Power consumption (VA)	UV Lamp	Purewater Tank
				W	D	H				
PX-0060 μ -X00	18.2 M Ω ·cm	\leq 5ppb	< 1.0L/min	290	428	598	Approx.24	60 (AC100-240V)	●	Built In
PX-0060 μ -000		—							—	

※The quality of treated water depends on the quality of raw water.

※TOC value is not displayed.